

Lens

OK

Save lenses during radiotherapy

- Efficient shielding
- Patient-comfortable
- High-quality CT images
- Safe silicone & replaceable



Lens OK

Save lenses during radiotherapy

Complete Therapeutic Lights
Uncomplicated Cancer Cure



Product Features

01. Tungsten Shield

Efficient shielding of electron beam

02. Transparent Silicone Cover

Easy to check the lens,
patient safety and convenience enhancement

03. Silicone of Sufficient Thickness

Reduce hot does due to
reduction of back-scattering

04. Replaceable Silicone

Long-term usable

05. Provide a Dedicated CT device

High-quality CT images

06. Fixed hole for Positioning


Prevent movement of shield during treatment

07. Eye-fitting Structure



Eye area treatment available(proton therapy, etc.)

Radiation Dose Distribution

| 6 MeV electron beam



| 9 MeV electron beam



RaphaRAD Inc.

165, Yeoksam-ro, Gangnam-gu, 06247 Seoul, Korea (Haesung Bldg)

TEL. +82-1811-7543

FAX. +82-2-567-7543

E-MAIL. info@rapharad.com

www.rapharad.com